



ETHOS X APPLICATION REPORTS

Microwave Green Extraction of Natural Products




MILESTONE
HELPING
CHEMISTS

Introduction

The microwave heating is generated by two mechanisms: ionic conduction and dipole rotation. The electromagnetic field applied by the microwave irradiation (MW) causes the migration of ions of the solvent and the resistance of this ion conduction results in friction that heats the solution. The electromagnetic field produces also a rearrangement of dipoles, which contributes to the heating up of the solvent as well. This conversion of electromagnetic energy into thermal energy is peculiar and unique for microwave heating process. In conventional processes the energy is transferred to the material by convection, conduction and radiation phenomena from the external material, while in MAE the microwave energy leads to directly molecular interactions with the electromagnetic field and consequent heating up of thesuscepting material^[1].

The microwave extraction (MAE) process occurs as the result of changes in the cell structure caused by electromagnetic waves. Its process acceleration and high extraction yield may be the result of a synergistic combination of two transport phenomena: heat and mass gradients working in the same direction. In that sense MAE differs from the conventional extraction methods (solid-liquid or simply extraction) where the mass transfer occurs from the inside to the outside, although the heat transfer occurs from the outside to the inside of the substrate.



Microwave selective heating

Furthermore in MAE the dissipation of the heating takes place volumetrically inside the irradiated medium, while in conventional extraction the heat is transferred from the heating medium to the interior of the sample. There are six main steps that characterize the microwave extraction process conducted with a solvent: 1) penetration of the solvent into the solid matrix; 2) solubilization and/or break down of components; 3) transport of the solute out of the solid matrix by internal diffusion; 4) migration of the extracted solute from the external surface of the solid into the bulk solution by external diffusion; 5) movement of the extract with respect to the solid; 6) separation and discharge of the extract and solid^[2]. Regarding the extraction process itself, three main subprocesses are known to take place and can be therefore pointed out: 1) desorption at an approximately constant velocity of the substrate from the outer surface of the particle (equilibrium phase); 2) after an intermediary transition phase, a mass transfer by internal diffusion and convection takes place in the solid-liquid interface; 3) in the last step the solute must overcome the interactions that bind it to the matrix and diffuse into the extracting solvent (external diffusion). This last step is considered as the limiting step of the process, since the extraction rate is really low and due to its irreversibility^[3].

Important parameters in Microwave-Assisted Extraction

EFFECT OF EXTRACTION TIME

Extraction time represents another parameter to be taken into consideration during MAE. The period of heating is indeed an important factor in MAE, since a selective heating occurs and extraction times are very short (from a few minutes to a half-hour) in order to avoid thermal degradation or oxidation of target compounds, that might be sensitive against overheating^[4]. The overheating depends on the dielectric properties of the solvent and it must be avoided especially in case of thermolabile constituents. In that sense, when longer extraction time is required, multiple steps extraction cycles are the best solution to enhance extraction yields, avoiding long heating^[4b, 5]. In flavonoids extraction from *R. astragali* for instance, it was found an increase in yield with time up to an exposure of 25 minutes and then the extraction yield started to decrease^[5].

EFFECT OF MICROWAVE POWER AND EXTRACTION TIME

Similarly to the period of heating, microwave power and consequently temperature are factors that must be properly set in order to maximize the extraction efficiency, without causing the degradation of thermally sensitive compounds. It is known that in MAE temperature is controlled by incident microwave power that controls the amount of energy provided to the matrix, which is converted to heat energy in the dielectric material^[2]. Microwave power is closely related to the extraction time required, the quantity of sample and of course to the type of sample we want to obtain the extract from. It acts as a driving force because it causes localized heating in plant matrix which allows the solute to diffuse out from the plant cells and to dissolve in the solvent. For that reason, an increase in the microwave power improves in general the extraction yield, ending up in shorter extraction time^[4b, 6]. On the other hand it must be however taken into consideration that a high microwave power might also cause a rapid breakage of the cell wall, ending up with the leaching out both of the target compounds and of undesired impurities, in addition to causing the degradation of thermally sensitive compounds^[7].

EFFECT OF THE SAMPLE CHARACTERISTICS

Another important parameter affecting the MAE process are the characteristics of the sample. It is well known that the higher the contact surface area is, the higher the extraction efficiency results, since finer particles allow a deeper penetration of the microwave^[8]. A prior grinding and homogenizing of the samples are therefore strongly recommended. In some cases soaking of the dried plant material in the extracting solvent before MAE has also resulted in improved yield^[7]. The moisture of the matrix acts as an extraction solvent as well, since it heats evaporates and causes internal ruptures of the cells, increasing though the extraction yield^[9].



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX00 - BASIL FRAGRANCES

Basil Solvent-Free Microwave
Extraction (SFME) and Microwave
Hydrodistillation (MWHD)



Introduction

Ocimum basilicum L. (Lamiaceae), respectively, named basil, is an aromatic herb that has been used traditionally as a medicinal herb in the treatment of headaches, coughs, diarrhea, constipation, warts, worms and kidney malfunctions. It has a long history as culinary herb, thanks to its foliage adding a distinctive flavor to many foods. It is also a source of aroma compounds and essential oils containing biologically active constituents that possess insecticidal, nematocidal, fungistatic and antimicrobial properties^[1].

[1] O. Politeo, M. Jukic, M. Milos, Food Chemistry 2007, 101, 379-385.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from basil can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHd) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHd techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHd techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 1*).

| Fresh Basil, SFME | | | | | | |
|-------------------|----------------------------------|-----------|---------|-------|------------------------|-----------|
| Reactor | Weighted fresh raw material [g] | Power [W] | Chiller | | Volatile fraction [mL] | Yield [%] |
| | | | 1kW | 2.1kW | | |
| Small | 500 | 500 | • | | 0.4 | 0.07 |
| Medium | 1580 | 1580 | | • | 1.3 | 0.08 |
| Large | 3720 | 1800 | | • | 3 | 0.08 |
| | | | | | | |
| Dry Basil, MWHD | | | | | | |
| Reactor | Weighted dry soaked material [g] | Power [W] | Chiller | | Volatile fraction [mL] | Yield [%] |
| | | | 1 kW | 2.1kW | | |
| Small | 500 | 500 | • | | 0.2 | 0.04 |
| Medium | 1580 | 1580 | | • | 0.5 | 0.03 |
| Large | 3720 | 1800 | | • | 1.5 | 0.04 |

Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Basil through SFME and MWHD. It is the unique modern concept of the

antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| No. | Identified compound | Peak area (%) | RI ^a HP-20M | RI ^a HP-101 |
|-----|-----------------------|---------------|------------------------|------------------------|
| 1 | β -Pinene | 0.1 | - | 949 |
| 2 | Limonene | 0.1 | 1180 | 1005 |
| 3 | 1,8-Cineole | 4.0 | 1185 | 1006 |
| 4 | Camphor | 0.5 | 1477 | 1109 |
| 5 | Linalool | 28.6 | 1518 | 1092 |
| 6 | Bornyl acetate | 0.5 | 1545 | 1252 |
| 7 | Terpinen-4-ol | 0.7 | 1563 | 1154 |
| 8 | α -Bergamotene | 2.2 | 1564 | 1407 |
| 9 | Caryophyllene | 0.3 | - | 1385 |
| 10 | Aloaromadendrene | 0.1 | - | 1450 |
| 11 | Estragole | 21.7 | 1632 | 1177 |
| 12 | α -Terpineol | 1.0 | 1653 | 1176 |
| 13 | Germacrene D | 0.3 | 1673 | 1444 |
| 14 | α -Humulene | 0.2 | - | 1417 |
| 15 | Carvone | 0.4 | 1685 | 1207 |
| 16 | β -Cubebene | 0.5 | 1694 | 1059 |
| 17 | β -Burbonene | t | - | 1354 |
| 18 | β -Elemene | 0.3 | - | 1364 |
| 19 | α -Cadinene | 0.2 | 1716 | 1426 |
| 20 | Calamenene | 0.2 | - | 1483 |
| 21 | α -Amorphene | 1.0 | 1710 | 1479 |
| 22 | β -Farnesene | 0.2 | - | 1452 |
| 23 | Δ -Cadinene | 0.1 | 1724 | 1486 |
| 24 | α -Bisabolene | 0.1 | - | 1506 |
| 25 | (Z)-Methyl cinnamate | 1.6 | 1900 | 1281 |

Table 2. Chemical composition of basil essential oil

| No. | Identified compound | Peak area (%) | RI ^a HP-20M | RI ^a HP-101 |
|-----|----------------------|---------------|------------------------|------------------------|
| 26 | Methyl eugenol | 3.1 | 1959 | 1378 |
| 27 | (E)-Methyl cinnamate | 14.3 | 2019 | 1364 |
| 28 | Spatulenol | 0.8 | 2066 | - |
| 29 | Eugenol | 5.9 | 2105 | 1368 |
| 30 | Carvacrol | t | 2118 | 1814 |
| 31 | α -Cadinol | 7.1 | 2120 | 1614 |
| 32 | Torreyol | 0.2 | 2173 | - |
| 33 | Chavicol | 0.7 | - | - ^b |
| | Total | 97.0 | | |

—, not identified.

t, trace (<0.1%).

^a Retention indices relative to C₈–C₂₂ n-alkanes on polar HP-20M and apolar HP-101 column.

^b Retention times is outside of retention times of homologous series of C₈–C₂₂ n-alkanes (identified by MS).

Table 2 (continued).




MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX02 - CITRUS PEEL FRAGRANCES

Citrus Peel Solvent-Free Microwave
Extraction (SFME) and Microwave
Hydrodistillation (MWHD)



Introduction

Citrus essential oils are the most widely used essential oils in the world. They are obtained as by-products of the citrus processing. They are used as aroma flavor in many food products, including alcoholic and non-alcoholic beverages, candy, gelatins. In pharmaceutical industries they are employed as flavoring agents to mask unpleasant tastes of drugs. In perfumery and cosmetic, they are used in many preparations. The traditional way to extract essential oils is by cold pressing the citrus peels. The oil is present in oil sacs or oil glands located at different depths in the peel and the cuticles of the fruit. Peel and cuticle oils are removed mechanically by cold pressing and since cold pressing yields a watery emulsion, this emulsion is then centrifuged to separate out the essential oil. Distillation is also used in some countries as an economical way to recover the oils. During distillation, the Citrus peels exposed to boiling water or steam, release their essential oils through evaporation. Researchers in many universities are working on novel techniques that could lead to compact, safe,

efficient, energy saving, and sustainable extraction processes. Solvent-Free Microwave Extraction (SFME) as upcoming extraction techniques have been reported for the extraction of fragrances and flavors from citrus peel^[3].

[3] N. Bousbia, M. Abert Vian, M. Ferhat, B. Meklati, F. Chemat, Journal of Food Engineering 2009, 90, 409-413.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from cistrus peel can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Citrus peel, SFME | | | | | |
|-------------------------|----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted fresh raw material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1kW | 2.1kW | |
| Small | 500 | 500 | • | | 1.0 |
| Medium | 1580 | 1580 | | • | 1.0 |
| Large | 3720 | 1800 | | • | 1.3 |

| Dry Citrus peel, MWHD | | | | | |
|-----------------------|-----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted dry soaked material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1 kW | 2.1kW | |
| Small | 500 | 500 | • | | 0.56 |
| Medium | 1580 | 1580 | | • | 0.57 |
| Large | 3720 | 1800 | | • | 0.61 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Citrus peel through SFME and MWH. It is the unique modern concept of

the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWH are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| No. | Compounds ^a | R.I. ^b | R.I. ^c | SFME |
|-----|--------------------------------|-------------------|-------------------|--------------|
| | Monoterpenes | | | 92.76 |
| 1 | Pinene<Alpha-> | 926 | 1023 | 1.75 |
| 2 | Pinene<Beta-> | 974 | 1109 | 15.35 |
| 3 | Myrcene<Beta-> | 988 | 1165 | 1.33 |
| 4 | Carene<Delta-3-> 1101 | 1101 | 1290 | 0.20 |
| 5 | Limonene | 1030 | 1206 | 65.25 |
| 6 | Terpinene<Gamma-> | 1103 | 1285 | 8.08 |
| | Oxygenated Monoterpenes | | | 92.76 |
| 7 | Linalool | 1125 | 1538 | 0.18 |
| 8 | Citronellal | 1167 | 1478 | 0.05 |
| 9 | Terpin-4-ol | 1191 | 1590 | 0.42 |
| 10 | Terpineol<Alpha-> | 1203 | 1677 | 0.56 |
| 11 | Nerol | 1237 | 1781 | 0.49 |
| 12 | Neral | 1268 | 1670 | 0.68 |
| 13 | Geraniol | 1271 | 1828 | 0.60 |
| 14 | Geranial | 1284 | 1714 | 0.89 |
| | Sesquiterpenes | | | 1.06 |
| 15 | Elemene<Beta-> | 1373 | 1583 | - |
| 16 | Caryophellene<E-> | 1391 | 1594 | 0.18 |
| 17 | Bergamotene<Alpha-Trans-> | 1437 | 1577 | 0.28 |
| 18 | Humulene<Alpha-> | 1450 | 1657 | 0.04 |
| 19 | Farnesene<(E)-Beta-> | 1453 | 1650 | 0.02 |
| 20 | Germacrene D | 1477 | 1696 | - |
| 21 | Valencene | 1488 | 1705 | 0.04 |
| 22 | Bisabolene<(Z)-Alpha-> | 1498 | 1761 | 0.03 |
| 23 | Bisabolene(Beta-) | 1508 | 1718 | 0.44 |

Table 3. Chemical compositions of essential oils from citrus peel obtained by SFME

| No. | Compounds ^a | R.I. ^b | R.I. ^c | SFME |
|-----|------------------------------------|-------------------|-------------------|-------------|
| | Oxygenated Sesquiterpenes | | | 0.03 |
| 24 | Elemol | 1540 | 1381 | - |
| 25 | Nerolidol<E-> | 1555 | 2026 | - |
| 25 | Bisabolol<Alpha-> | 1684 | 2212 | 0.03 |
| 27 | Nootkatone | 1799 | 2250 | - |
| | Other oxygenated compounds | | | 0.86 |
| 28 | Nonanal<N-> | 1126 | 1400 | 0.06 |
| 29 | Citronellyl Acetate | 1342 | 1645 | 0.04 |
| 30 | Neryl Acetate | 1351 | 1706 | 0.19 |
| | Extraction time (min) | | | 180 |
| | Yield (%) | | | 0.8 |
| | Total oxygenated compounds (%) | | | 4.78 |
| | Total non oxygenated compounds (%) | | | 93.82 |

^a Essential oil compounds sorted by chemical families and percentages calculated by GC-FID on non-polar HP5MS™ capillary column.

^b Retention indices calculated on non-polar HP5MS™ capillary column.

^c Retention indices calculated on polar Carbowax™-PEG capillary column.

Table 3 (continued).



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX03 - FRANKINCENSE FRAGRANCES

Hydrodistillation (MWHD)



Introduction

Gum resins from *Boswellia* species, also known as frankincense, have been used as a major ingredient in Ayurvedic and Chinese medicine to treat a variety of health-related conditions. Both frankincense chemical extracts and essential oil prepared from *Boswellia* species gum resins exhibit anti-neoplastic activity, and have been investigated as potential anti-cancer agents. The anti-cancer activity is mediated through multiple signaling pathways. In addition, frankincense essential oil overcomes multicellular resistant and invasive phenotypes of human breast cancer cells. Fast and green Extraction of Frankincense essential oil turns out therefore to be extremely important. This essential oil is obtained from a resin from the bark of a shrub originally from the area surrounding the Red Sea, in Somalia and Arabia. To collect the resin, fine incisions are made in the bark, and drops of sap appear and dry in large, odorous yellow droplets^[1].

[1] X. Ni, M. Suhail, Q. Yang, A. Cao, K.-M. Fung, R. Postier, C. Woolley, G. Young, J. Zhang, H.-K. Lin, BMC Complementary and Alternative Medicine 2012, 12.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Frankincense can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) technique allows the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 1*).

| Frankincense, MWHD | | | | | |
|--------------------|---|-----------|---------|-------|------------------------|
| Reactor | Weighted dry material + added water [g]* | Power [W] | Chiller | | Volatile fraction [mL] |
| | | | 1 kW | 2.1kW | |
| Small | 1050 | 1050 | • | | 5.5 |
| Medium | 3318 | 1800 | | • | 17.6 |
| Large | 7816 | 1800 | | • | 43.8 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Frankincense resin through SFME and MWHD. It is the unique modern

concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| RRI | Compounds | Boswellia rivae (%) | | | |
|------|-------------------------------|---------------------|------|----------------------------------|------|
| 1032 | α -Pinene | 5.3 | 1597 | Bornyl acetate | 0.4 |
| 1035 | α -Thujene | 1.3 | 1600 | β -Elemene | 0.3 |
| 1076 | Camphene | 0.1 | 1611 | Terpinen-4-ol | 1.4 |
| 1118 | β -Pinene | 0.6 | 1524 | trans-Dihydrocarvone | 0.2 |
| 1132 | Sabinene | 1.2 | 1639 | trans-p-Mentha-2,8-diene-1-ol | 3.9 |
| 1159 | δ -3-Carene | 9.6 | 1642 | Thuj-3-en-10-al | 0.2 |
| 1187 | o-Cymene | 2.5 | 1648 | Myrtenal | 0.5 |
| 1203 | Limonene | 14.8 | 1651 | Sabinaketone | 0.2 |
| 1213 | 1,8-Cineole | 0.3 | 1657 | Umbellulone | 0.1 |
| 1266 | (E)- β -Ocimene | 0.4 | 1663 | cis-Verbenol | 0.5 |
| 1278 | m-Cymene | 0.4 | 1664 | trans-Pinocarveol | 2.2 |
| 1280 | p-Cymene | 2.9 | 1678 | cis-p-Mentha-2,8-diene-1-ol | 0.9 |
| 1424 | o-Methylanisol | 0.2 | 1683 | trans-Verbenol | 6.8 |
| 1430 | α -Thujone | 0.1 | 1700 | p-Mentha-1,8-diene-4-ol | 0.4 |
| 1439 | γ -Campholene aldehyde | 0.1 | 1706 | α -Terpineol | 1.4 |
| 1444 | 2,5-Dimethylstyrene | 0.2 | 1709 | α -Terpinyl acetate | 1.0 |
| 1450 | trans-Linalool oxide | 0.1 | 1720 | trans-Sabinol | 0.3 |
| 1451 | β -Thujone | 0.7 | 1725 | Verbenone | 4.3 |
| 1458 | cis-1,2-Limonene epoxide | 4.6 | 1751 | Carvone | 1.6 |
| 1468 | trans-1,2-Limonene epoxide | 0.5 | 1804 | Myrtenol | 0.7 |
| 1474 | trans-Sabinene hydrate | 0.9 | 1811 | trans-p-Mentha-1(7),8-diene-2-ol | 0.3 |
| 1478 | cis-Linalool oxide | 0.1 | 1845 | trans-Carveol | 2.5 |
| 1498 | (E)- β -Ocimene epoxide | 0.2 | 1856 | m-Cymen-8-ol | 3.1 |
| 1499 | α -Campholene aldehyde | 0.9 | 1864 | p-Cymen-8-ol | 2.0 |
| 1536 | Pinocamphone | 0.3 | 1882 | cis-Carveol | 0.7 |
| 1553 | Linalool | 0.2 | 1896 | cis-p-Mentha-1(7),8-diene-2-ol | 0.2 |
| 1556 | cis-Sabinene hydrate | 0.9 | 1949 | Piperitenone | 0.6 |
| 1565 | 8,9-Limonene epoxide-I | 0.6 | 2113 | Cumin alcohol | 0.1 |
| 1571 | 8,9-Limonene epoxide-II | 0.6 | 2198 | Thymol | 0.1 |
| 1586 | Pinocarvone | 0.5 | 2239 | Carvacrol | 0.1 |
| | | | | Total | 88.1 |

Table 2. The composition of Boswellia Rivae essential oil



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX05 - GARLIC FRAGRANCES

Garlic Solvent-Free Microwave
Extraction (SFME) and Microwave
Hydrodistillation (MWHD)



Introduction

Garlic (*Allium sativum*) belongs to a group of dietary supplements that lessen the incidence of cardiovascular and cerebrovascular diseases by reducing cholesterol concentration. The beneficial effect of garlic on health confirm several studies which were showed that garlic has been evaluated for lowering blood pressure, cholesterol and glucose concentration, reduce blood lipids as well as for the prevention of arteriosclerosis and cancer. The biological activities of garlic including antibacterial, parasiticial, antithrombotic, antioxidant and antidiabetic actions have been known for a long time. The unique flavor and health-promoting functions of garlic are generally attributed to its rich content of sulfur-containing compounds: alliin, g-glutamylcysteine and their derivatives.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from garlic can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Garlic, SFME | | | | | |
|--------------------|---------------------------------|-----------|---------|-------|------------------------|
| Reactor | Weighted fresh raw material [g] | Power [W] | Chiller | | Volatile fraction [mL] |
| | | | 1kW | 2.1kW | |
| Small | 1000 | 1000 | • | | 2 |
| Medium | 3160 | 1800 | | • | 8.2 |
| Large | 7445 | 1800 | | • | 22.3 |

Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Garlic through SFME and MWHD. It is the unique modern concept of the

antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| Compounds | LRI _{HP1} | LRI _{INNO} | SFME (% ± SD) | Identification methods |
|---|--------------------|---------------------|---------------|------------------------|
| Dimethyl sulfide | - | 750 | tr | SM, LRI |
| Dimethyl disulfide | 731 | 1099 | tr | SM, LRI, Std |
| Methyl ethyl disulfide | 818 | 1119 | 0.5 | SM, LRI |
| Diallyl sulfide | 840 | 1150 | 0.5 | SM, LRI, Std |
| Methyl allyl disulfide | 894 | 1290 | 0.8 ± 0.1 | SM, LRI, Std |
| (Z)-prop-1-enyl methyl disulfide | 917 | 1251 | tr | SM, LRI |
| (E)-prop-1-enyl methyl disulfide | 919 | 1275 | 0.1 | SM, LRI |
| Dimethyl trisulfide | 941 | 1340 | 0.2 | SM, LRI, Std |
| Diallyl disulfide | 1056 | 1491 | 26.3 ± 0.9 | SM, LRI, Std |
| Allyl (Z)-prop-1-enyl disulfide | 1073 | 1390 | 2.9 ± 0.1 | SM, LRI |
| Allyl (E)-prop-1-enyl disulfide | 1082 | 1415 | 7.7 ± 0.1 | SM, LRI |
| Allyl methyl trisulfide | 1115 | 1601 | 7.9 ± 0.1 | SM, LRI, Std |
| Methyl (E)-propenyl trisulfide | 1138 | - | 0.1 | Tentative |
| 3-vinyl-(4H)-1,2-dithiin | 1156 | 1735 | 1.1 | SM, LRI |
| Unknown 1 | 1169 | - | 0.1 | - |
| 2-vinyl-(4H)-1,3-dithiin | 1178 | 1830 | 3.1 | SM, LRI |
| Diallyle trisulfide | 1285 | 1825 | 39.7 ± 0.8 | SM, LRI, Std |
| Propyl propenyl trisulfide ^d | 1290 | 1781 | 0.3 | - |
| Allyl propenyl trisulfide | 1300 | 1798 | 2.0 | SM, LRI |
| 3,5-diethyl-1,2,4-trithiolane | 1321 | 1788 | 0.2 | SM, LRI |
| Unknown 2 | 1348 | 2011 | 0.4 | - |
| Diallyl tetrasulfide | 1508 | - | 1.5 | SM, LRI, Std |
| 2,4-dimethyl-5,6-dithia-2,7-nonadienal | 1718 | > 2400 | 0.9 ± 0.1 | SM, LRI |
| 2,4-dimethyl-5,6-dithia-2,7-nonadienal | 1730 | > 2400 | 0.4 ± 0.1 | SM, LRI |
| Unknown 3 | 1775 | - | 0.1 | - |

^a Compounds are listed in order of their elution time from a HP-1 column. Compositional values less than 0.1% are denoted as traces (tr). Presence of a compound is indicated by its GC-FID percentage with S.D., absence is indicated by "-".

^b RI = retention indices are determined on HP-1 and INNOWAX column using the homologous series of n-alkanes (C6-C24).

^c S.D. = standard deviation.

^d Correct isomer not identified

tentative: tentatively identified by MS and RI without standard compound co-injection

Unknown 1: 186 (M+, 7.1); 162(23.2); 100 (15); 97 (28.3); 73(50.1); 60 (16.0); 56(15.6); 59 (28.2); 57 (31.8); 45 (45.5); 41(100).

Unknown 2: 178 (M+, 3.4); 172(16.6); 170 (100); 128(42.3); 106(25.3); 73(12.3); 64(62.3); 59(17.3); 45(26.8); 42 (11.7); 41(31).

Unknown 3: 179(M+, 5.6); 147(89.0); 106(10.5); 105(33.4); 75(14.6); 73(98.7); 64 (14.8); 57(10.7); 47 (15.4); 45(29.7); 41(100).

Table 3. Essential oil from garlic (*Allium sativum* oil) obtained by SFME



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX07 - GINGER FRAGRANCES

Ginger Solvent-Free Microwave
Extraction (SFME) and Microwave
Hydrodistillation (MWHD)



Introduction

Due to its composition in valuable natural compounds, ginger represents a suitable matrix for extraction of essential oil as well as many other bioactive compounds. It contains indeed products of interest such as essential oils (1–4%), phenolics (gingerols and 6-shogaol, 1–2%), and total carbohydrates (60–75%). Ginger, and more specifically rhizomes are variously used as food product or in traditional medicine. In the food industry, rhizomes are mainly used for spices or condiments (fresh or dried), candy or as juice after cold mechanical pressing. Due to the fact that mechanical pressing does not alter the chemical composition of the pressed product, this process provides huge amounts of press cake still containing high amounts of bioactive compounds, which is currently considered as waste. The novelty of the extraction of these products through MHG and SFME relies on the extraction of compounds achieved without addition of solvent or water. The only water used in the process is the matrix water naturally present in the plant cells ^[1].

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Ginger can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Ginger, SFME | | | | | |
|--------------------|---------------------------------|-----------|---------|-------|------------------------|
| Reactor | Weighted fresh raw material [g] | Power [W] | Chiller | | Volatile fraction [mL] |
| | | | 1kW | 2.1kW | |
| Small | 1000 | 1000 | • | | 0.7 |
| Medium | 3160 | 1800 | | • | 2.2 |
| Large | 7445 | 1800 | | • | 5.2 |

Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Ginger through SFME and MWHD. It is the unique modern concept of the

antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| | GR | GP | 0.6 W g ⁻¹ | 0.8 W g ⁻¹ | 1.0 W g ⁻¹ | 1.2 W g ⁻¹ | 1.4 W g ⁻¹ | 1.6 W g ⁻¹ | 1.8 W g ⁻¹ |
|---|------|------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Essential oil Yield (g per 100 g fresh plant material). Major compounds (%) | 0.2 | 0.3 | 0.2 | 0.2 | 0.2 | 0.2 | 0.2 | 0.2 | 0.2 |
| β-Pinene | 1.2 | 1.0 | 2.3 | 2.6 | 2.4 | 2.6 | 2.3 | 2.4 | 2.2 |
| Camphene | 4.3 | 3.8 | 9.1 | 10.3 | 9.2 | 10.0 | 9.1 | 9.4 | 9.1 |
| Sabinene | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 | 0.1 |
| Sulcatone | 0.0 | 0.8 | 1.2 | 2.8 | 3.3 | 3.2 | 3.0 | 3.2 | 2.9 |
| Myrcene | 0.6 | 0.6 | 0.0 | 1.4 | 1.4 | 1.4 | 1.3 | 1.3 | 1.1 |
| α-Phellandrene | 0.2 | 0.1 | 0.2 | 0.3 | 0.3 | 0.3 | 0.3 | 0.3 | 0.1 |
| Limonene | 0.9 | 0.9 | 1.7 | 1.9 | 1.9 | 1.9 | 1.7 | 1.8 | 1.7 |
| β-Phellandrene | 4.6 | 4.2 | 8.7 | 10.4 | 10.3 | 10.2 | 9.7 | 10.0 | 8.6 |
| Terpinolene | 0.1 | 0.1 | 0.1 | 0.1 | 0.3 | 0.3 | 0.3 | 0.3 | 0.2 |
| Linalol | 0.2 | 0.2 | 0.3 | 0.4 | 0.4 | 0.4 | 0.4 | 0.4 | 0.4 |
| Borneol | 0.5 | 0.6 | 0.8 | 0.9 | 1.0 | 0.9 | 1.0 | 1.0 | 1.1 |
| α-Terpineol | 0.2 | 0.3 | 0.4 | 0.5 | 0.5 | 0.5 | 0.5 | 0.5 | 0.6 |
| Citronellol | 0.1 | 0.3 | 0.2 | 0.5 | 0.4 | 0.3 | 0.4 | 0.4 | 0.8 |
| Neral | 1.7 | 0.5 | 0.4 | 1.3 | 1.5 | 1.7 | 1.5 | 1.5 | 1.3 |
| Geraniol | 0.1 | 0.2 | 0.1 | 0.3 | 0.3 | 0.2 | 0.2 | 0.2 | 0.6 |
| Geranial | 3.3 | 1.0 | 0.6 | 1.9 | 2.2 | 2.6 | 2.3 | 2.5 | 2.3 |
| Geranyl acetate | 0.3 | 0.1 | 0.4 | 0.2 | 0.2 | 0.2 | 0.2 | 0.2 | 0.2 |
| α-Curcumene | 3.5 | 13.9 | 17.0 | 7.6 | 7.2 | 6.6 | 7.0 | 6.8 | 9.9 |
| Germacrene D | 1.6 | 1.3 | 0.1 | 1.3 | 1.4 | 1.4 | 1.4 | 1.4 | 0.7 |
| Zingiberene | 35.7 | 25.2 | 18.4 | 23.2 | 24.0 | 24.0 | 25.1 | 24.3 | 18.4 |
| α-Farnesene | 6.5 | 6.5 | 6.3 | 5.4 | 5.5 | 5.5 | 5.7 | 5.5 | 5.7 |
| β-Bisabolene | 5.7 | 6.8 | 0.0 | 4.8 | 4.7 | 4.6 | 4.8 | 4.7 | 5.4 |
| β-Sesquiphellandrene | 12.1 | 13.9 | 12.3 | 9.9 | 9.9 | 9.7 | 10.2 | 9.8 | 10.4 |
| Antioxidants Total content (g per 100 g plant material DW). Major compounds (g per 100 g plant material DW) | 1.17 | 0.90 | 0.57 | 1.24 | 1.06 | 1.18 | 1.22 | 1.37 | 1.18 |
| 6-Gingerol | 0.77 | 0.58 | 0.31 | 0.81 | 0.65 | 0.79 | 0.81 | 0.92 | 0.79 |
| 8-Gingerol | 0.15 | 0.11 | 0.07 | 0.14 | 0.11 | 0.14 | 0.14 | 0.17 | 0.14 |
| 10-Gingerol | 0.23 | 0.19 | 0.11 | 0.18 | 0.19 | 0.19 | 0.19 | 0.21 | 0.19 |
| 6-Shogaol | 0.02 | 0.02 | 0.08 | 0.11 | 0.10 | 0.08 | 0.09 | 0.08 | 0.08 |

Table 3. Volatile compounds and antioxidants extracted from ginger plant material by SFME



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX09 - LAVANDER FRAGRANCES

Lavender Solvent-Free Microwave Extraction (SFME) and Microwave Hydrodistillation (MWHD)



Introduction

Lavender is one of the most useful medicinal plants. Commercially, the lavender provides several important essential oils to the fragrance industry, including soaps, colognes, perfumes, skin lotions and other cosmetics. In food manufacturing, lavender essential oil is employed in flavouring beverages, ice cream, candy, baked goods and chewing gum. The essential oils of *Lavandula* species are obtained by steam distillation of the fresh flowering spikes. Oil quality is assessed by oil chemical composition and by the organoleptic opinion. In addition, a large range of medical uses for this plant have also been reported. These include antispasmodic, sedative, antihypertensive, antiseptic, healing and anti-inflammatory properties, all of which render it highly appreciated in phytotherapy and aromatherapy. In food manufacturing, lavender essential oil has been employed in flavoring beverages, ice cream, baked goods and chewing gum ^[1].

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Lavander can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Lavander, SFME | | | | | | |
|----------------------|-----------------------------------|-----------|---------|-------|------------------------|-----------|
| Reactor | Weighted fresh raw material [g]* | Power [W] | Chiller | | Volatile fraction [mL] | Yield [%] |
| | | | 1kW | 2.1kW | | |
| Small | 1000 | 1000 | • | | 2.5 | 0.25 |
| Medium | 3160 | 1800 | | • | 6.7 | 0.21 |
| Large | 7445 | 1800 | | • | 15 | 0.22 |
| | | | | | | |
| Dry Lavander, MWHD | | | | | | |
| Reactor | Weighted dry soaked material [g]* | Power [W] | Chiller | | Volatile fraction [mL] | Yield [%] |
| | | | 1 kW | 2.1kW | | |
| Small | 1000 | 1000 | • | | 1.7 | 0.17 |
| Medium | 3160 | 1800 | | • | 6.6 | 0.21 |
| Large | 7445 | 1800 | | • | 14.2 | 0.19 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Lavander through SFME and MWHD. It is the unique modern concept of the

antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| No. | Compounds | R.R.I. ^a | MWHD ^b (%) |
|-----|---------------------------|---------------------|-----------------------|
| | Monoterpenes | | 3.54 |
| 1 | α -Thujene | 907 | 0.08 |
| 2 | α -Pinene | 914 | 0.51 |
| 3 | Camphene | 933 | 0.32 |
| 4 | Sabinene | 963 | 0.14 |
| 5 | β -Pinene | 968 | 0.59 |
| 6 | β -Myrcene | 985 | 0.50 |
| 7 | 3-Carene | 1008 | 0.22 |
| 8 | Limonene | 1024 | tr. |
| 9 | (Z)- β -Ocimene | 1031 | 0.33 |
| 10 | (E)- β -Ocimene | 1040 | 0.37 |
| 11 | γ -Terpinene | 1050 | 0.09 |
| 12 | Terpinolene | 1078 | 0.37 |
| | Oxygenated monoterpenes | | 78.29 |
| 13 | 1,8-Cineole | 1027 | 7.23 |
| 14 | Sabinene hydrate-cis | 1058 | 0.66 |
| 15 | Linalool | 1099 | 47.82 |
| 16 | Camphor | 1137 | 11.82 |
| 17 | Borneol | 1161 | 4.15 |
| 18 | Terpin-4-ol | 1174 | 5.94 |
| 19 | p -Cymen-8-ol | 1179 | tr. |
| 20 | α -Terpineol | 1186 | 0.68 |
| | Sesquiterpenes | | 2.77 |
| 21 | α -Bergamotene cis | 1400 | 0.10 |
| 22 | β -Caryophyllene | 1412 | 1.28 |

Table 3. Yields, extraction times, grouped compounds and chemical compositions of essential oils obtained by MWHD from lavender flowers

| No. | Compounds | R.R.I. ^a | MWHD ^b (%) |
|-----|--------------------------------------|---------------------|-----------------------|
| 23 | α -Santalene | 1414 | 0.15 |
| 24 | (E)- β -Farnesene | 1453 | 0.63 |
| 25 | Sesquiterpene 1 | 1474 | 0.61 |
| | Oxygenated sesquiterpenes | | 0.29 |
| 26 | Caryophyllene oxide | 1573 | 0.11 |
| 27 | α -Bisabolol | 1677 | 0.18 |
| | Other oxygenated compounds | | 15.01 |
| 28 | Octan-3-one | 977 | 0.78 |
| 29 | Octan-3-ol | 990 | 0.26 |
| 30 | Dihydromyrcenol | 1063 | 0.34 |
| 31 | n.i. | 1141 | 0.37 |
| 32 | n.i. | 1188 | 2.00 |
| 33 | n.i. | 1232 | 0.43 |
| 34 | Linalool acetate | 1254 | 10.74 |
| 35 | Geranyl acetate | 1377 | 0.08 |
| | Yield (%) | | 8.86 |
| | Total extraction time (min) | | 10 |
| | Heating time from 20 to 100 °C (min) | | 5 |
| | Real extraction time (min) | | 5 |

tr., trace; n.i., non-identified.

^a R.R.I., relative retention indices relative to C₈–C₂₂ n-alkanes on SBP5™ capillary column.

^b MWHD, microwave hydrodistillation.

Table 3 (continued).



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX11 - MINT FRAGRANCES

Mint Solvent-Free Microwave
Extraction (SFME) and Microwave
Hydrodistillation (MWHD)



Introduction

The mint species have a great importance, both medicinal and commercial. Indeed, leaves, flowers and stems of *Mentha* spp. are frequently used in herbal teas or as additives in commercial spice mixtures for many foods to offer aroma and flavor. In addition, *Mentha* spp. has been used as a folk remedy for treatment of nausea, bronchitis, flatulence, anorexia, ulcerative colitis, and liver complaints due to its antiinflammatory, carminative, antiemetic, diaphoretic, antispasmodic, analgesic, stimulant, emmenagogue, and anticatharrhal activities. Commercially, the most important mint species are peppermint (*M. x piperita*), spearmint (*M. spicata*), and corn mint (*M. canadensis*). From these species, corn mint is cultivated only because of oil production. Peppermint oil is one of the most popular and widely used essential oils, mostly because of its main components menthol and menthone. Corn mint is the richest source of natural menthol (Sharma and Tyagi 1991; Shasany et al. 2000). Carvone-scented mint plants, such as spearmint, are rich in carvone and are widely used as spices and

cultivated in several countries. Peppermint oil is used for flavouring pharmaceuticals and oral preparations, such as toothpastes, dental creams, and mouth washes. It is also used as a flavouring agent in cough drops, chewing gums, confectionery and alcoholic liqueurs. It is used in medicines for internal use. Its pleasant taste makes it an excellent gastric stimulant ^[1].

[1] H. Hajlaoui, N. Trabelsi, E. Noumi, M. Snoussi, H. Fallah, R. Ksouri, A. Bakhrouf, World Journal of Microbiology and Biotechnology 2009, 25, 2227-2238.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Mint can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Mint, SFME | | | | | | |
|------------------|-----------------------------------|-----------|---------|-------|------------------------|-----------|
| Reactor | Weighted fresh raw material [g]* | Power [W] | Chiller | | Volatile fraction [mL] | Yield [%] |
| | | | 1kW | 2.1kW | | |
| Small | 500 | 500 | • | | 4.8 | 0.98 |
| Medium | 1580 | 1580 | | • | 16 | 1.01 |
| Large | 3720 | 1800 | | • | 37.2 | 1.0 |
| | | | | | | |
| Dry Mint, MWHD | | | | | | |
| Reactor | Weighted dry soaked material [g]* | Power [W] | Chiller | | Volatile fraction [mL] | Yield [%] |
| | | | 1 kW | 2.1kW | | |
| Small | 500 | 500 | • | | 4.3 | 0.86 |
| Medium | 1580 | 1580 | | • | 13.7 | 0.87 |
| Large | 3720 | 1800 | | • | 30.9 | 0.83 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Mint through SFME and MWHD. It is the unique modern concept of the antiquated

Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large) complying with a high range of extraction-scale needs.

| | Fresh leaves SFME (g Kg ⁻¹) | Dried leaves MWHD (g Kg ⁻¹) |
|-----------------------|---|---|
| α -Pinene | 4.0 | 7.4 |
| Sabinene | 4.6 | 6.3 |
| β -Pinene | 7.2 | 10.8 |
| 2-Thujene | 15.9 | 14.7 |
| 3-Octanol | - | 2.6 |
| Limonene | 67.5 | 8.67 |
| 1,8-Cineole | 24.8 | 3.45 |
| (E)- β -Ocimene | 6.3 | 5.2 |
| (Z)- β -Ocimene | 3.4 | 2.0 |
| γ -Terpinene | 3.6 | 1.1 |
| 3-Carene | 8.1 | 11.1 |
| p-Menth-1-en-8-ol | - | 1.4 |
| (-)-4-Terpineol | 7.1 | 2.6 |
| t-Dihydrocarvone | 18.9 | 17.3 |
| c-Carveol | 7.1 | 3.3 |
| t-Carveol | - | 2.3 |
| D-Carvone | 602.7 | 601.7 |
| Piperitone | 4.9 | 8.2 |
| t-Carvone oxide | 2.8 | 4.0 |
| c-Carvone oxide | - | 4.2 |
| Dihydroedulan II | - | 1.0 |
| Dihydroedulan I | - | 1.5 |
| Isolimonene | 5.8 | - |
| Dihydrocarvyl acetate | - | 3.2 |
| t-Carveyl acetate | - | 1.8 |
| β -Bourbonene | 29.7 | 24.7 |
| β -Elemene | 4.2 | 5.5 |
| Isocaryophyllene | - | 1.1 |

Table 3. Main components of *M. spicata* L. var. *rubra* EO yield (g kg⁻¹)

| | Fresh leaves SFME (g Kg ⁻¹) | Dried leaves MWHD (g Kg ⁻¹) |
|---------------------------------|---|---|
| β-Caryophyllene | 44.8 | 42.3 |
| α-Caryophyllene | 3.9 | 3.8 |
| (+)-Epi-bicyclosquiphellandrene | 9.8 | 9.0 |
| Germacrene D | 17.3 | 16.1 |
| Bicyclogermacrene | 4.4 | 3.5 |
| α-Murolene | - | 1.3 |
| Calamenene | 9.3 | 5.6 |
| Caryophyllene oxide | 16.8 | 3.0 |
| 1,4-Cadinadiene | 5.9 | 2.2 |
| τ-Murolol | 4.1 | 1.2 |
| Total | 945.0 | 954.0 |

Table 3 (continued).



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX13 - ORANGE PEEL FRAGRANCES

Orange Peel Solvent-Free
Microwave Extraction (SFME) and
Microwave Hydrodistillation (MWHD)



Introduction

Orange essential oil is used to confer the aroma and flavor of orange to a wide variety of products, such as carbonated drinks, ice-creams, cakes, air-fresheners, perfumes and so on. They are also being used in the design of new products, to which they add aroma and flavor 2, 3 and 4. Another application takes advantage of the germicide properties of some of their components. In that sense, a small amount of d-limonene was very effective in the germicide treatment of waste waters. Carotenoid pigments, found in orange extracts, are important for health, not only because of their nutritional value as precursors of vitamin A, but also because of their antioxidant potential. They also seem to have certain anticarcinogenic properties. These compounds are used in food coloring [1].

[1] B. Mira, M. Blasco, A. Berna, S. Subirats, Journal of Supercritical Fluids 1999, 14, 95-104.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Orange peel can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHd) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHd techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHd techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Orange peel, SFME | | | | | |
|-------------------------|----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted fresh raw material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1kW | 2.1kW | |
| Small | 500 | 500 | • | | 4.9 |
| Medium | 1580 | 1580 | | • | 16 |
| Large | 3720 | 1800 | | • | 69 |

| Dry Orange peel, MWHd | | | | | |
|-----------------------|-----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted dry soaked material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1 kW | 2.1kW | |
| Small | 500 | 500 | • | | 2.8 |
| Medium | 1580 | 1580 | | • | 9 |
| Large | 3720 | 1800 | | • | 23.1 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Orange peel through SFME and MWHD. It is the unique modern concept of

the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| N° | Compound | I ^a | SFME (%) |
|-------------------------|-----------------------|----------------|----------|
| Monoterpenes | | | |
| 1 | α -Pinene | 928 | 0.60 |
| 2 | Sabinene | 968 | 0.23 |
| 3 | β -Myrcene | 988 | 1.81 |
| 4 | 3-Carene | 1001 | tr |
| 5 | δ -3-Carene | 1007 | 0.05 |
| 6 | Limonene | 1034 | 96.20 |
| 7 | α -Terpinolene | 1087 | 0.01 |
| Oxygenated monoterpenes | | | |
| 8 | Linalool | 1093 | 0.17 |
| 9 | Trans-Limonene oxide | 1135 | 0.02 |
| 10 | Citronellal | 1150 | 0.02 |
| 11 | β -Citronellol | 1227 | 0.06 |
| Sesquiterpene | | | |
| 12 | α -Copaene | 1372 | 0.01 |
| 13 | β -Cubebene | 1382 | tr |
| 14 | β -elemene | 1386 | 0.01 |
| 15 | Caryophyllene (E) | 1417 | 0.01 |
| 16 | α -Humulene | 1452 | 0.05 |
| 17 | Germacrene-D | 1479 | tr |
| 18 | Valencene | 1490 | 0.01 |
| 19 | Germacrene-A | 1503 | 0.01 |
| 20 | δ -Cadinene | 1520 | 0.01 |

Table 3. Chemical composition of orange EOs obtained by SFME

| N° | Compound | I ^a | SFME (%) |
|----------------------------|----------------------|----------------|----------|
| Oxygenated sesquiterpenes | | | |
| 21 | Caryophyllene oxide | 1589 | tr |
| 22 | Cis, trans -Farnesol | 1694 | - |
| 23 | α -Sinensal | 1754 | tr |
| Other oxygenated compounds | | | |
| 24 | Decanal | 1203 | 0.16 |
| 25 | n-Dodecanal | 1404 | 0.02 |
| Extraction time (min.) | | 6 | 120 |
| Yield (%) | | 5.43 | 5.45 |

tr. tracesb0.01%.

^a I, Retention indices relative to C5–C28 n-alkanes calculated on non-polar HP5MS capillary column. Percentages calculated by GC–FID on non-polar HP5MS capillary column. Essential oil compounds sorted by chemical families.

Table 3 (continued).



MILESTONE

H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX14 - ROSE FRAGRANCES

Rosa Damascena Mill. Solvent-Free
Microwave Extraction (SFME) and
Microwave Hydrodistillation (MWHd)



Introduction

Roses have been used since the earliest times in rituals, cosmetics, perfumes, medicines and aromatherapy. A great variety of garden roses also exist, which are bred less for fragrance and more for color and shape. Even with the high price of roses and the advent of organic synthesis, rose oils are still the most widely used essential oils in perfumery. In fact, because of the labor-intensive production process and the low content of oil in rose blooms, rose oil is very expensive and is often called as 'liquid gold'. For the production of essential rose oil as well as rose extracts, the two rose species most often used are *Rosa damascena* (pink damask rose) and *Rosa centifolia* (light pink cottage rose). Whereas the former is predominantly used for rose oil production, the oil of the latter is usually extracted with solvents such as petroleum ether or n-hexane in order to obtain rose concrete. Rose concrete, the result of solvent extraction, is mainly composed of fragrance-related substances, and contains large quantities of paraffins, fatty acids, fatty acid methyl esters, di- and tri-terpenic

compounds and pigments. The concrete can be therefore processed through Solvent-Free Microwave Extraction (SFME) as upcoming extraction technique for the extraction of essential oil to eliminate non-volatile compounds such as paraffins, shorten extraction time, reduce organic solvent consumption, improve extraction yield, enhance quality of the extract, prevent pollution and reduce sample preparation costs ^[1].

[1] M. Mohamadi, T. Shamspur, A. Mostafavi, Journal of Essential Oil Research 2012, 25, 55-61.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Rosa Damascena Mill. can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave

Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques. (See Cookbook for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 1*).

| Fresh Rosa Damascena Mill., SFME | | | | | |
|----------------------------------|---------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted fresh raw material [g] | Power [W] | Chiller | | Yield [%] |
| | | | 1kW | 2.1kW | |
| Small | 400 | 400 | • | | 0.1 |
| Medium | 1264 | 1264 | | • | 0.14 |
| Large | 2978 | 1800 | | • | 0.13 |

Time, Power
 The extractions were carried out till complete recovery of the fragrance
 ≤ 1800 g: Power(W) = Weight(g).
 > 1800g: Power = 1800W
 Chiller settings:
 ≤ 900g, 1 kW Chiller
 > 900g, 2.1 kW Chiller

Table 1

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Dry Rosa Damascena Mill. through SFME and MWHD. It is the unique modern

concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| Compound | RI ^a | SFME (% ^b) |
|--------------------------------|-----------------|------------------------|
| Monoterpenes | | 0.6 |
| α -Pinene | 858 | 0.2 |
| Sabinene | 896 | tr ^c |
| β -Pinene | 898 | 0.1 |
| β -Myrcene | 916 | 0.1 |
| α -Terpinene | 937 | tr |
| Limonene | 949 | tr |
| γ -Terpinene | 977 | tr |
| α -Terpinolene | 1003 | tr |
| Oxygenated monoterpenes | | 49.3 |
| Rosefuran | 1009 | tr |
| Linalool | 1017 | 0.4 |
| <i>trans</i> -Rose oxide | 1023 | 0.3 |
| Citronellal | 1056 | tr |
| Neroloxide | 1061 | tr |
| 4-Terpineol | 1081 | 0.1 |
| β -Citronellol | 1157 | 34.7 |
| Geraniol | 1173 | 9.1 |
| Eugenol | 1262 | 1.5 |
| Geranyl acetate | 1596 | 0.6 |
| Methyl eugenol | 1309 | 2.3 |
| Neryl acetate | 1926 | 0.2 |
| Sesquiterpenes | | 4.1 |
| β -Bourbonene | 1302 | 0.2 |
| <i>trans</i> -caryophyllene | 1326 | 0.5 |

Table 2. Qualitative and quantitative composition of rose essential oils obtained by SFME

| Compound | RI ^a | SFME (% ^b) |
|-------------------------------------|-----------------|------------------------|
| α -Guaiene | 1340 | 0.4 |
| α -Humulene | 1348 | 0.5 |
| Germacrene D | 1368 | 1.4 |
| β -Selinene | 1369 | 0.1 |
| Caryophyllen (1I) | 1378 | 0.1 |
| α -Selinene | 1379 | - |
| δ -Guaiene | 1381 | 0.4 |
| (<i>E,E</i>)- α -Farnesene | 1385 | 0.1 |
| δ -Cadinene | 1394 | tr |
| Ledene | 1518 | 0.3 |
| Oxygenated sesquiterpenes | | 2.2 |
| Elmol | 1419 | tr |
| Nerolidol | 1437 | 0.1 |
| Caryophyllene oxide | 1446 | tr |
| γ -Eudesmol | 1505 | 0.1 |
| β -Eudesmol | 1516 | 0.2 |
| <i>cis</i> -Farnesol | 1555 | 1.8 |
| Hydrocarbons | | 31.4 |
| Pentadecane | 1388 | 0.2 |
| Hexadecane | 1492 | 0.1 |
| 8-Heptadecene | 1547 | 0.2 |
| Heptadecane | 1569 | 1.8 |
| <i>cis</i> -9-Tricosene | 1630 | 0.2 |
| Octadecane | 1669 | 0.4 |
| (<i>Z</i>)-5-Nonadecene | 1737 | 3.5 |
| Nonadecane | 1765 | 15.1 |
| (<i>E</i>)-9-Eicosene | 1811 | 0.6 |
| Eicosane | 1836 | 1.2 |
| 9-Nonadecene | 1882 | 0.2 |
| 1-Nonadecene | 1894 | 0.2 |
| <i>n</i> -Heneicosan | 1910 | 6.0 |
| Docosan | 1972 | 0.1 |
| <i>cis</i> -9-Tricosene | 2029 | 0.2 |
| Tricosane | 2039 | 0.9 |
| Tetracosane | 2099 | tr |
| Pentacosane | 2160 | 0.2 |
| <i>n</i> -Heptacosane | 2274 | 0.2 |
| Other oxygenated compounds | | 0.6 |
| Phenylethyl alcohol | 813 | 0.4 |
| Heptanal | 1011 | tr |

Table 2 (continued).

| Compound | RI ^a | SFME (% ^b) |
|---------------------------------------|-----------------|------------------------|
| Tetradecanal | 1024 | 0.1 |
| Nonanal | 1485 | tr |
| Total oxygenated compounds | | 52.1 |
| Total non-oxygenated compounds | | 36.1 |
| Extraction time (min) | | 210 |
| Oil yield (w/w%) | | 0.018 |

^aRetention indices relative to C6–C27 *n*-alkanes on HP-1MS column.

^b%, relative percentage obtained on HP-1 column using GC/MS detector.

^ctr, 60.05. *Significant at $p \leq 0.05$ based on *F*-value determined by analysis of variance (ANOVA). **Significant at $p \leq 0.01$ based on *F*-value determined by ANOVA. CV, coefficient of variation.

Table 2 (continued).



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX16 - ROSEMARY FRAGRANCES

Rosemary Solvent-Free Microwave Extraction (SFME) and Microwave Hydrodistillation (MWHD)



Introduction

Rosemary (*Rosmarinus officinalis* L.) is an important herb on the world food and aromatherapy market. Natural antioxidants such as those present in rosemary essential oil may be an alternative source for compounds capable of protecting lipids in foods. The essential oil secreted by glandular trichomes is mainly located in leaves. Rosemary essential oil is also used as an antibacterial and antifungal agent. Nevertheless, it has been noted that these activities often depend on the origin of the rosemary plant and the method of extraction. Since both of these quality parameters can greatly influence the chemical composition of rosemary oil, Solvent-Free Microwave Extraction (SFME) as upcoming extraction techniques have been reported for the extraction of fragrances from Rosemary ^[2].

[2] N. Tigrine-Kordjani, B. Meklati, F. Chemat, International Journal of Aromatherapy 2006, 16, 141-147.

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Rosemary can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

| Fresh Rosemary, SFME | | | | | |
|----------------------|----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted fresh raw material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1kW | 2.1kW | |
| Small | 500 | 500 | • | | 0.63 |
| Medium | 1580 | 1580 | | • | 0.65 |
| Large | 3720 | 1800 | | • | 0.72 |

| Dry Rosemary, MWHD | | | | | |
|--------------------|-----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted dry soaked material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1 kW | 2.1kW | |
| Small | 500 | 500 | • | | 0.38 |
| Medium | 1580 | 1580 | | • | 0.38 |
| Large | 3720 | 1800 | | • | 0.4 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Rosemary through SFME and MWHF. It is the unique modern concept of the

antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHF are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

| N° | Compound ^a | SFME (%) | RI ^b | RI ^c |
|----|--------------------------|----------|-----------------|-----------------|
| | Monoterpene hydrocarbons | 68.56 | | |
| 1 | Tricyclene | 0.26 | 921 | 1011 |
| 2 | α -Pinene | 44.05 | 936 | 1023 |
| 3 | Camphene | 6.14 | 951 | 1103 |
| 4 | Verbenene | 0.77 | 955 | 1121 |
| 5 | β -Pinene | 2.61 | 980 | 1109 |
| 6 | Myrcene | 1.94 | 995 | 1149 |
| 7 | α -Phellandrene | 0.31 | 995 | 1165 |
| 8 | γ -3-Carene | 0.08 | 1014 | 1290 |
| 9 | α -Terpinene | 0.86 | 1020 | 1083 |
| 10 | para-Cymene | 1.27 | 1025 | 1250 |
| 11 | Limonene | 5.48 | 1030 | 1206 |
| 12 | γ -Terpinene | 3.08 | 1052 | 1251 |
| 13 | Terpinolene | 1.71 | 1092 | 1287 |
| | Oxygenated monoterpenes | 24.87 | | |
| 14 | Linalool | 2.00 | 1106 | 1538 |
| 15 | α -Campholenal | 1.24 | 1122 | 1471 |
| 16 | Camphor | 7.82 | 1149 | 1514 |
| 17 | Pinocarvone | 1.33 | 1160 | 1548 |
| 18 | Borneol | 2.57 | 1173 | 1679 |
| 19 | Terpin-4-ol | 2.07 | 1184 | 1590 |
| 20 | α -Terpineol | 0.77 | 1198 | 1677 |
| 21 | Verbenone | 6.37 | 1207 | 1696 |
| 22 | Geraniol | 0.70 | 1279 | 1828 |

Table 3. Chemical composition of *Rosmarinus officinalis* essential oils obtained by SFME.

| N° | Compound ^a | SFME (%) | RI ^b | RI ^c |
|----|--------------------------------|-----------|-----------------|-----------------|
| | Sesquiterpene hydrocarbons | 1.91 | | |
| 23 | E-Caryophyllene | 0.95 | 1425 | 1470 |
| 24 | α -Humulene | 0.42 | 1450 | 1657 |
| 25 | γ -Curcumene | 0.04 | 1469 | 1738 |
| 26 | β -Bisabolene | 0.43 | 1508 | 1714 |
| 27 | β -Sesquiphellandrene | 0.07 | 1519 | 1776 |
| | Oxygenated sesquiterpenes | 0.26 | | |
| 28 | Caryophyllene oxide | 0.10 | 1570 | 1977 |
| 29 | α -Bisabolol | 0.16 | 1684 | 2022 |
| | Other oxygenated compounds | 1.03 | | |
| 30 | Bornyl acetate | 0.81 | 1263 | 1579 |
| 31 | Methyl eugenol | 0.12 | 1397 | 2032 |
| 32 | Z-Methyl jasmonate | 0.10 | 1635 | 2349 |
| | | | | |
| | Extraction time (min) | 180 | | |
| | Yield (%) | 0.35±0.07 | | |
| | Total oxygenated compounds | 26.16 | | |
| | Total non-oxygenated compounds | 70.47 | | |

^a Essential oil compounds sorted by chemical families and percentages calculated by GC–FID on non-polar HP5MSTM capillary column.

^b Retention indices relative to C₅–C₂₈ n-alkanes calculated on non-polar HP5MSTM capillary column.

^c Retention indices relative to C₅–C₂₈ n-alkanes calculated on polar CarbowaxTM-PEG capillary column.

Table 3 (continued).



MILESTONE
H E L P I N G
C H E M I S T S

APPLICATION REPORT

EX17 - THYMUS FRAGRANCES

Thymus vulgaris Solvent-Free
Microwave Extraction (SFME) and
Microwave Hydrodistillation (MWHD)



Introduction

Thymus vulgaris L. (common thyme), a member of the Labiatae family, is an aromatic/medicinal plant of increasing economic importance for North America, Europe, North Africa and Asia. Thyme is one of many aromatic plants that has been utilized in variety of food products to provide a flavor specific to this herb. Studies indicating the antiseptic, carminative, antimicrobial, and antioxidative properties of thyme have also been published. From the medicinal point of view, thyme has been used as a culinary herb and also as a herbal medicine. The essential oils of thyme are responsible for the typical spicy aroma of the plant leaves. These oils are stored in glandular peltate trichomes situated on both sides of the leaves. Results published on the chemical composition of thyme oil revealed that most of the oil was produced from flowering plants. In the plant's life cycle, the oil production is usually at its highest level during this period. For *T. vulgaris* and *T. pulegioides*, such finding was reported in the early 1960s ^[1].

Why to choose Microwave Fragrances set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from *Thymus vulgaris* can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the SFME or MWHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

Results and experimental procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 1*).

| Fresh <i>Thymus vulgaris</i> , SFME | | | | | |
|-------------------------------------|----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted fresh raw material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1kW | 2.1kW | |
| Small | 400 | 400 | • | | 0.32 |
| Medium | 1264 | 1264 | | • | 0.35 |
| Large | 2978 | 1800 | | • | 0.37 |

| Dry <i>Thymus vulgaris</i> , MWHD | | | | | |
|-----------------------------------|-----------------------------------|-----------|---------|-------|-----------|
| Reactor | Weighted dry soaked material [g]* | Power [W] | Chiller | | Yield [%] |
| | | | 1 kW | 2.1kW | |
| Small | 400 | 400 | • | | 0.16 |
| Medium | 1264 | 1264 | | • | 0.15 |
| Large | 2978 | 1800 | | • | 0.18 |

*Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

Important remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please take care to seal properly the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from *Thymus vulgaris* through SFME and MWHD. It is the unique modern concept of

the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large) complying with a high range of extraction-scale needs.

| N° | RT ^b (min) | Compound | I _k ^c | Relative peak area ^a (%) SFME |
|----|-----------------------|------------------------|-----------------------------|--|
| 1 | 16.05 | α -Thujene | 930 | 0.53 \pm 0.02 |
| 2 | 16.53 | α -Pinene | 938 | 0.86 \pm 0.01 |
| 3 | 17.31 | Camphene | 952 | 0.53 \pm 0.01 |
| 4 | 19.04 | 1-Octen-3-ol | 983 | 2.64 \pm 0.31 |
| 5 | 19.60 | β -Myrcene | 993 | 1.30 \pm 0.17 |
| 6 | 19.91 | 3-Octanol | 999 | 0.19 \pm 0.03 |
| 7 | 20.39 | α -Phellandrene | 1008 | 0.18 \pm 0.01 |
| 8 | 20.71 | Δ -3-Carene | 1014 | 0.09 \pm 0.01 |
| 9 | 21.12 | α -Terpinene | 1021 | 1.73 \pm 0.14 |
| 10 | 21.73 | p -Cymene | 1033 | 17.57 \pm 0.78 |
| 11 | 21.99 | 1,8-Cineole | 1037 | 1.31 \pm 0.12 |
| 12 | 23.55 | γ -Terpinene | 1066 | 8.54 \pm 0.02 |
| 13 | 23.91 | Trans-Sabinene hydrate | 1073 | 0.94 \pm 0.05 |
| 14 | 24.95 | Terpinolene | 1093 | 0.27 \pm 0.05 |
| 15 | 25.59 | Linalool | 1105 | 2.43 \pm 0.27 |
| 16 | 29.20 | Borneol | 1176 | 1.11 \pm 0.21 |
| 17 | 29.69 | Endo-Borneol | 1185 | 1.41 \pm 0.21 |
| 18 | 29.80 | Terpinen-4-ol | 1187 | 0.63 \pm 0.16 |
| 19 | 31.20 | α -Terpineol | 1216 | 0.17 \pm 0.00 |
| 20 | 32.34 | Methyl thymylether | 1240 | 0.14 \pm 0.07 |
| 21 | 34.32 | Geraniol | 1281 | 0.39 \pm 0.07 |
| 22 | 35.90 | Thymol | 1315 | 40.20 \pm 3.03 |
| 23 | 36.50 | Carvacrol | 1328 | 6.84 \pm 0.68 |

Table 2. Chemical compositions of essential oils obtained by solvent-free microwave extraction (SFME) of thyme aerial parts

| N° | RT ^b (min) | Compound | I _K ^c | Relative peak area ^a (%) HD |
|-----------------------------|-----------------------|---------------------|-----------------------------|--|
| 24 | 38.22 | Thymyl acetate | 1366 | 0.16 ± 0.03 |
| 25 | 41.36 | β-Caryophyllene | 1438 | 2.86 ± 0.27 |
| 26 | 42.72 | α-Humulene | 1470 | 0.64 ± 0.22 |
| 27 | 42.98 | Geranyl acetate | 1477 | 0.42 ± 0.09 |
| 28 | 45.44 | Δ-Cadinene | 1536 | 0.40 ± 0.05 |
| 29 | 48.17 | Caryophyllene oxide | 1604 | 1.42 ± 0.21 |
| Total peak area (%) | | | 94.31 | 95.91 |
| Total extraction time (min) | | | 120 | 240 |
| Yield (%) | | | 2.39 ± 0.06% | 2.52 ± 0.00% |

^a Mean ± SD (n = 2).

^b Retention time.

^c Kovats Retention Index (I_K) relative to C₉–C₁₈ n-alkanes on the HP-5MS column.

Table 2 (continued).